

The Genetics of Athletic Performance

DNA Sport

Sasha Mannion, MSc

Daniel Meyersfeld, PhD

Ian Craig, MSc, CSCS, INLPTA

2017



Table of Contents

1. Genetic variants associated with injury susceptibility	2
1.1. Collagen 1 Alpha 1 (COL1A1): Sp1 G>T	2
1.2. Collagen 5 Alpha 1 (COL5A1): BstUI C>T	3
1.3. Growth Differentiation Factor 5 (GDF5): 5' UTR C>T	4
2. Genetic variants associated with recovery from exercise	5
2.1. Interleukin-6 (IL-6): -174 G>C.....	5
2.2. Interleukin 6 Receptor (IL-6R): Asp358Ala (A>C)	6
2.3. C-Reactive Protein (CRP): 219 G>A	7
2.4. Tumour Necrosis Factor Alpha (TNFA): -308 G>A	8
2.5. Superoxide Dismutase 2 (SOD2): 208 T>C (Val16Ala)	8
2.6. Endothelial Nitric Oxide Synthase (eNOS): Asp298Glu (G>T)	9
3. Genetic variants associated with athletic performance	10
3.1. Angiotensinogen (AGT): Met238Thr (T>C).....	10
3.2. Angiotensin I-Converting Enzyme (ACE): Insertion (I)/Deletion (D).....	11
3.3. Bradykinin Receptor B2 (BDKRB2): Exon 1 +9/-9 (C>T).....	12
3.4. Vascular Endothelial Growth Factor (VEGF): -634 C>G	13
3.5. Nuclear Respiratory Factor 2 (NRF-2); GABPB1: A>G.....	13
3.6. Peroxisome Proliferator-Activated Receptor-Gamma Coactivator-1 (PPARGC1A): Gly482Ser (G>A)	14
3.7. Peroxisome Proliferator-Activated Receptor-Alpha (PPARA): Intron 7 G>C	15
3.8. Beta 2 Adrenergic Receptor (ADRB2): Arg16Gly (A>G); Gln27Glu (C>G)	15
3.9. Thyrotropin Releasing Hormone Receptor (TRHR): rs7832552 C>T	16
3.10. Alpha-actinin 3 (ACTN3): 577 R>X.....	17
3.11. Vitamin D Receptor (VDR): Taq1 T>C	18
3.12. Cytochrome P450 1A2 (CYP1A2): -163 C>A	18
3.13. Circadian Locomotor Output Cycles Kaput (CLOCK): 3111 T>C	19
4. References	20

1. Genetic variants associated with injury susceptibility

Soft tissue injuries are a common occurrence in competitive and recreational sports and activities. These injuries can be due to spontaneous ruptures (acute onset injuries) or can be chronic overuse injuries due to repetitive strain. Injury to any soft tissue is considered a multifactorial condition that is caused by a complex interaction between intrinsic risk factors, including genetic risk factors, and extrinsic, non-genetic, risk factors.

By the identification of genetic risk factors and the various biological processes that contribute to injury development, therapeutic strategies for injury prevention and treatment can be established.

1.1. Collagen 1 Alpha 1 (COL1A1): Sp1 G>T

Type I collagen is the most abundant protein in the body and is a major structural component of bones, tendons, ligaments and joint cartilage. Type I collagen is a heterodimer that consists of two $\alpha 1$ chains and one $\alpha 2$ chain, encoded by the *COL1A1* and *COL1A2* genes respectively. Along with small amounts of type III and type V collagen, type I collagen molecules assemble into a highly organised array that results in a mature connective tissue structure, essential for tensile strength.¹

Variations within the *COL1A1* gene have been linked to risk of soft tissue injury as well as alterations in bone mineral density, osteoarthritis and fracture risk.

The rare TT genotype of the functional *COL1A1* Sp1 G>T SNP has been linked to a substantially reduced risk of cruciate ligament and shoulder dislocation injuries in comparison to the GG genotype.² In agreement, Posthumus *et al.* (2009) showed that the TT genotype was significantly under-represented in an injured group (with anterior cruciate ligament ruptures) in comparison to matched active controls,³ demonstrating the possible protective role of this genotype. Combined analysis of the *COL1A1* Sp1 G>T SNP from independent studies further advocates the protective role of the TT genotype in acute soft tissue ruptures,²⁻⁵ indicating the importance of this SNP being included in a model to determine injury risk.

Functional analysis shows that when the *COL1A1* Sp1 T allele is present, an increase in the affinity for the transcription factor Sp1 results in increased *COL1A1* gene expression.⁶ The consequence of this is not yet fully understood, however it is hypothesized that an increased production of COL1A1 in TT individuals increases the tensile strength of tendons and ligaments, while the resulting imbalance between the $\alpha 1$ and $\alpha 2$ chains conversely impairs bone strength.^{2,6}

1.2. Collagen 5 Alpha 1 (COL5A1): BstUI C>T

COL5A1 encodes the $\alpha 1$ chain of type V collagen, a minor constituent of soft tissues such as tendons and ligaments that is essential for life. Type V collagen is a fibrillary collagen that intercalates with type I collagen to form heterotypic fibres. Along with other proteins, type V collagen plays an important structural role in regulating fibrillogenesis (fibril formation) and fibril diameter within tendons, ligaments and other connective tissues.⁷

The *COL5A1* BstUI C>T polymorphism has been associated with injuries such as Achilles tendinopathy and anterior cruciate ligament ruptures as well as range of motion and endurance performance.⁸

The CC genotype of the *COL5A1* BstUI polymorphism has been associated with protection from soft tissue injuries. In 2006 Mokone *et al.* found that the CC genotype was associated with protection from chronic Achilles tendinopathy, with these individuals being less likely to develop the symptoms associated with the injury.⁹ This study was replicated in South African and Australian populations yielding the same results.¹⁰ When this polymorphism was investigated as a risk factor for anterior cruciate ligament ruptures, the CC genotype was significantly over-represented in active female controls when compared to the ACL injured counterparts.¹¹ More recently the *COL5A1* BstUI T allele has been linked to an increased risk of tennis elbow.¹² In 2017, Brown *et al.* further verified the role of *COL5A1* as a genetic determinant of injury risk.¹³

The function of the *COL5A1* polymorphism has been associated with altered mRNA stability, affecting *COL5A1* gene expression. The T allele is linked to increased *COL5A1* mRNA stability and the C allele to lower synthesis of the alpha 1 chain of type V collagen.¹⁴ It is suggested that the relative content of type V collagen in soft tissues alters the mechanical properties and injury susceptibility.⁸ Animal models have reported an association with the increase in content of type V collagen and a decrease in fibril diameter and biomechanical properties of tendons.¹⁵

COL5A1 BstUI C>T has also been associated with range of motion.¹⁶ A 2010 study by Brown *et al.* showed that older individuals with the CC genotype had a greater range of motion and protection from injury.¹⁷ Due to the association between increased stiffness (decreased range of motion) and improved running economy, the *COL5A1* gene was investigated for an association with athletic performance. *COL5A1* Sp1 TT individuals were found to have a faster marathon completion in an Ironman triathlon in comparison to CT and CC individuals.¹⁸ An additional study showed that TT individuals were less flexible and completed a 56km ultramarathon significantly faster than participants with TC and CC genotypes.¹⁹

From the collective studies, the hypothesis is that the TT genotype results increased production of *COL5A1*, altered mechanical properties of musculoskeletal soft tissues, increased risk of injury, reduced joint ROM (flexibility), and increased endurance running ability.⁸

1.3. Growth Differentiation Factor 5 (GDF5): 5' UTR C>T

GDF5 encodes the Growth Differentiation Factor 5 protein, part of the bone morphogenetic protein (BMP) family and transforming growth factor β (TGF- β) superfamily. GDF5 is involved in tissue growth, development, homeostasis and healing, including that of bones, cartilage and various musculoskeletal soft tissues.²⁰ This gene has been linked to mechanical strength, with animal models showing significantly weaker tendons with lowered collagen content in GDF5 deficient mice.²⁰ Conversely, the addition of exogenous GDF5 showed improved ligament and tendon repair, tensile strength and biomechanical testing in rat and canine models.^{21,22}

The functional polymorphism within the 5'UTR of GDF5 has been shown to influence susceptibility to soft tissue injuries as well other multifactorial disorders.

In 2010 Posthumus *et al.* investigated the association of GDF5 with risk of developing Achilles tendinopathy. Individuals with the TT genotypes were identified as having a 2-fold increased risk of injury.²³ The T allele of this GDF5 polymorphism is correlated with reduced expression of the GDF5 gene within soft tissues. In line with this, the collective studies suggest that the reduced expression of this gene with the GDF5 T allele leads to an increased risk of exercise related soft tissue injuries.²⁴

The GDF5 T allele has also been strongly linked to an increased risk of osteoarthritis, critically of the knee.²⁴⁻²⁷

2. Genetic variants associated with recovery from exercise

The recovery genes included in DNA Sport are the best example within this genetic panel of the need to integrate training and nutritional advice when supporting an athlete's health and fitness. Two important biological areas are considered with regards to recovery: inflammation and oxidative stress.

Although not generally considered within sports science or sports nutrition, systemic inflammation and oxidative stress that results from heavy training and other lifestyle stresses, can greatly affect an athlete's health and rate of recovery between training bouts and sessions. Knowledge of the genes that code for these variables can provide athletes with an enormous opportunity to improve their specific health and training responses via nutrition, stress management and nutraceuticals. Six genetic polymorphisms for recovery are included in the DNA Sport test: they have been identified as those with strong scientific credibility and those which could make a difference to the training, nutrition and lifestyle interventions of an athlete, based on knowledge of their genes.

2.1. Interleukin-6 (IL-6): -174 G>C

IL-6 is a cytokine secreted by T cells and macrophages, which stimulates an immune response to trauma such as strenuous exercise, leading to inflammation in the muscle and fatty tissue. IL-6 stimulates energy mobilisation, causing an increase in body temperature. IL-6 is also a myokine (muscle cytokine) and regulates the gene encoding C-reactive protein (CRP). The immune system of athletes is affected by the intensity and duration of exercise training - the inevitable muscle damage results in an inflammatory repair process, which is mediated by inflammatory cytokines such as IL-6. Overtraining syndrome (OTS) has been hypothesised to be caused by excess cytokine release during exercise, resulting in a chronic inflammatory state.

The C allele has been linked to an increased IL-6 and CRP inflammatory response to exercise, which may induce more pronounced fatigue and prolong recovery times. Perhaps because IL-6 levels are hard to measure in the plasma, there was some discrepancy in the literature about whether the C or G allele is associated with higher IL-6 levels, but the majority and best evidence points towards the C allele.

In a study of 54 military recruits during their 8-week basic training, it was noted that acute exercise increased plasma IL-6 levels more in subjects with the CG genotype compared to GG homozygotes.²⁸ Interestingly, subjects with the CG genotype also made the greatest gains in VO₂max. In support of this, Yamin *et al.* (2008) found that subjects with the C allele experienced greater creatine kinase (CK) levels than GG homozygotes (3-times the risk of a massive CK response) after eccentric training, showing that the C allele may be a risk factor for exercise-induced muscle injury.²⁹ The IL-6 gene has been shown to be amongst a number of genes that were predictive of sports performance,³⁰ and the G allele seems to be the most common in power athletes compared to control and endurance athletes,³¹ which may be due to improved muscle repair after eccentric exercise.²⁹

Under-performance syndrome (UPS), another name for overtraining syndrome, is thought to be influenced by cytokine sickness, an over-production and/or intolerance to IL-6 and other cytokines.³²

This suggestion was supported by Robson-Ansley *et al.* (2007) who noted that an acute period of intensified training can suppress the innate immune system and chronically increase IL-6 levels.³³ These elevated cytokines can in-turn increase fatigue and malaise, which are related to the cytokine theories of UPS. The IL-6 genotype can also influence glucose levels: McKenzie *et al* (2004) found that baseline fasting glucose levels were higher in the CC genotype compared to carriers of the G allele: 6 months of aerobic training successfully decreased the glucose area under the curve during an oral glucose tolerance test, but a significant decrease only occurred in the GG genotype.³⁴

2.2. Interleukin 6 Receptor (IL-6R): Asp358Ala (A>C)

Interleukin 6 receptor (IL-6R) is a type I cytokine receptor, made up of a protein complex consisting of an IL-6 receptor subunit (IL-6R) and a signal transducer (gp130). The IL6R gene specifically encodes this IL-6R subunit, which in-turn influences IL-6 cytokine action. Proteolytic cleavage of the membrane-bound IL-6R protein leads to the generation of a soluble form of IL-6R (sIL-6R), which is able to bind to IL-6. The resulting IL-6/sIL-6R complex is also capable of binding to gp130 and inducing intracellular signalling,³⁵ increasing the biological activity and half-life of IL-6.³⁶ The resulting IL-6-induced inflammatory response has been linked to fatigue during exercise, the ability to recover from training sessions and potentially the risk of overtraining. Regular exercise reduces baseline long-term inflammatory biomarkers, but single high-intensity exercise sessions will acutely increase local inflammation.

An A to C SNP on the IL-6R gene has been shown to have an effect on sIL-6R and consequently IL-6 levels. Subjects with the CC or AC genotypes may have higher levels of sIL-6R and IL-6, thereby increasing the acute inflammatory effects of exercise. Individuals with these genotypes are generally advised to increase their recovery time between training bouts and training sessions, regularly check their inflammatory biomarkers and increase anti-inflammatory nutrition support.

It has been demonstrated in a study of 70 subjects that C allele carriers for the IL-6R gene have significantly higher baseline sIL-6R levels.³⁷ In African- and European-American subjects, one copy of the C allele has been shown to cause a 1.06 to 1.15-fold increase, and two copies of the C allele caused a 1.2 to 1.43-fold increase in IL-6 levels.³⁸ It has been clearly shown that IL-6 levels are increased during acute exercise bouts: during one hour of cycling at 90% of lactate threshold, Gray *et al* (2008, 2009) found a 5-fold increase in IL-6 and 1.2-fold increase in sIL-6R immediately afterwards. These levels returned to baseline 1.5 hours after exercise.^{36,39} Robson-Ansley *et al* (2011) also showed that the plasma concentrations of IL-6 increased post-exercise when athletes ran at 60% of their vVO₂max for 2 hours, followed by a 5km time trial. In this study, when a carbohydrate drink was ingested during the exercise, post-exercise IL-6 levels were significantly reduced.³³ The same investigators tracked 13 cyclists who rode 468 km over a six-day period: IL-6 levels were increased post-exercise on day 1, but remained unchanged for the following five days. sIL-6R, CRP and creatine kinase were unchanged post-exercise on day 1, but were elevated at baseline for the rest of the trial. It was demonstrated that sIL-6R was significantly correlated with CRP and that these inflammatory markers may affect subjective sensations of post-exercise fatigue.³³ The investigators concluded that strenuous,

prolonged exercise stimulated an acute-phase inflammatory response, which was maintained throughout the 6-day event.

2.3. C-Reactive Protein (CRP): 219 G>A

CRP is an acute-phase protein, which rises in response to inflammation in the body: IL-6 and other inflammatory cytokines trigger the synthesis of CRP and fibrinogen by the liver. Its physiological role is to bind to phosphocholine, which is expressed on the surface of dead or dying cells in order to activate the complement system, which is part of the innate immune system. CRP is used as a biomarker for acute or chronic inflammation. Regular exercise and favourable dietary habits can reduce baseline long-term inflammatory biomarkers (such as CRP), whereas single high-intensity exercise sessions will acutely increase local inflammation.

The CRP AA genotype, which is only found in 12% of the population, has been associated with lower levels of CRP compared to carriers of the G allele, who make up the majority of the population (88%).

Obisesan *et al* (2004) have shown that AA homozygotes have significantly lower CRP levels (40%) than G allele carriers and that CRP levels decrease significantly with 24 weeks of exercise training and a low-fat diet (although GG genotypes still had the highest CRP levels).⁴⁰ This has been confirmed by Eiriksdottir *et al* (2009), who found that carriers of the G allele had significantly higher CRP levels than non-carriers.⁴¹ Kullo *et al* (2007) have shown that higher circulating levels of IL-6, CRP and fibronectin are associated with lower VO₂max levels after controlling for other variables.⁴² This observation was confirmed by Kuo *et al* (2007) in a huge study of 1438 healthy adults.⁴³ When monitoring the recovery from an Ironman triathlon event, Neubauer *et al* (2008) found that CRP levels were raised post-exercise, along with several other immune and inflammatory markers. By 19 days post-event, most markers were back to normal except for CRP and myoglobin, which were still slightly raised, indicating low-grade inflammation for several days after such a large exercise bout.⁴⁴

It has been shown by a number of studies that the inflammation induced by eccentric exercise is greater during recovery when subjects eat a high-carbohydrate compared to a low-carbohydrate diet.^{45,46} For example, Depner *et al* (2010) took 12 subjects through 6 sets of 10 high-intensity eccentric contractions and found that post-exercise inflammation was reduced on a low-carbohydrate diet,⁴⁵ although it is important to remember that a carbohydrate drink during exhaustive exercise has been linked with a decrease in IL-6 levels.⁴⁷

There has been much research examining the use of antioxidants to decrease inflammation. For example, Phillips *et al* (2003) found that after eccentric exercise, subjects consuming antioxidant and fatty acid supplements had significantly lower levels of IL-6 and CRP.⁴⁸ However free radicals are involved in signalling roles during training adaptations and so antioxidant supplements may weaken these signals. There is currently much controversy in the antioxidant literature, so caution is generally recommended around high-dose supplementation long-term.⁴⁹

2.4. Tumour Necrosis Factor Alpha (TNFA): -308 G>A

TNF codes for the pro-inflammatory cytokine, tumor necrosis factor (TNF). It is produced chiefly by activated macrophages and is a member of a group of cytokines that stimulate the acute phase inflammatory reaction. The primary role of TNF is in the regulation of immune cells: it is able to induce fever, apoptotic cell death, sepsis (through IL-1 & IL-6 production), inflammation and to inhibit tumour growth and viral replication. Like other inflammatory cytokines, TNF levels increase during and after intensive exercise.

A, the minor allele on the TNFA gene, has been shown to cause 2-fold greater levels of transcription compared to the G form. The A allele is therefore associated with increased circulating TNF levels. The AA genotype, although very rare (2%), has been associated with increases in CRP levels during a 20-week exercise programme: CRP levels are generally expected to decrease with regular exercise training. An A allele carrier, with greater TNF levels, may experience increased levels of fatigue and poorer recovery times.

Like other inflammatory cytokines, TNF has been shown to be increased during and after highly-intensity exercise. During 3 hours of cycling or inclined walking, plasma concentrations of IL-1 β , IL-6 and TNF peaked at the end of exercise.⁵⁰ IL-1 β and TNF were shown to be still elevated 24 hours later. Lakka *et al* (2006) noted that prior to a 20-week exercise programme, the AA genotype had greater baseline CRP levels. What's more, the exercise programme decreased CRP levels less in the AA individuals than other genotypes.⁵¹ Nicklas *et al* (2005) took 213 older or overweight individuals with knee osteoarthritis through an 18 month walking and weight training exercise trial. It was found that walking distance and stair-climb time were better for individuals homozygous for the G allele.⁵² Nieman *et al.* (2005) demonstrated that carbohydrate intake during 2.5 hours of cycling at 60% max power decreased plasma cortisol, adrenaline, IL-6 and IL-10 responses, but did not affect TNF alpha levels.⁵³

2.5. Superoxide Dismutase 2 (SOD2): 208 T>C (Val16Ala)

The SOD2 gene codes for manganese superoxide dismutase (SOD2), which is a potent free radical scavenger within the cell, especially the mitochondria. The SOD2 enzyme converts superoxide free radicals to hydrogen peroxide. The mitochondria is commonly referred to as the workhorse of the cell and is therefore the site of many oxidative reactions during energy production, so free radicals are generated here. There is evidence that oxidation within the cell contributes to muscular fatigue and extreme exercise can cause increased lipid peroxidation and depleted levels of vitamin E. Long-term training can increase base levels of SOD2, whereas short-term bursts of intense activity will increase oxidative stress.

The SOD2 C allele has been associated with high oxidative stress biomarkers and people with this genotype who consume low levels of fruit and vegetables are at increased risk of developing long-

term disease. These levels of oxidative stress can be reduced by a diet that contains antioxidants from fruit and vegetables and supporting SOD2 activity.

In a study of 231 healthy students, it was found that SOD2 enzyme activity was 33% higher in CT or TT individuals compared to CC individuals and on average, SOD2 activity was 15% higher in females than males.⁵⁴ It has been shown that exercise training favourably increases baseline levels of antioxidant enzymes. Garcia-Lopez *et al.* (2007) demonstrated that mRNA levels of catalase, glutathione peroxidase (GP), MnSOD (mitochondrial) and CuZnSOD2 (cytosolic) were increased after 21 weeks of strength training, while endurance training only increased SOD2 and GP mRNA levels.⁵⁵ Mastaloudis *et al* (2004a) studied 22 runners during a 50 km ultramarathon. There were two groups: placebo or antioxidant (1g Vit C, 300mg Vit E).⁵⁶ F2-isoprostanes, a measure of lipid peroxidation, which were at similar levels between groups at baseline, increased during the run only in the placebo group. In the same group,⁵⁷ DNA damage increased mid-race, but returned to baseline 2-hours post-race.

2.6. Endothelial Nitric Oxide Synthase (eNOS): Asp298Glu (G>T)

The endothelium-derived nitric oxide (NO) plays a key role in the regulation of vascular tone and peripheral vascular resistance. It also has vasoprotective effects by suppressing platelet aggregation, leukocyte adhesion and smooth muscle cell proliferation.⁵⁸ "Impaired endothelial function, either as a consequence of reduced production/release or increased inactivation of endothelium-derived vasodilators, as well as interactions of NO with angiotensin, reactive oxygen species and oxidized lipoproteins, has detrimental functional consequences and is one of the most important cardiovascular risk factors".⁵⁸

The T allele has been associated with decreased activity of the eNOS enzyme which is linked to an increase in the presence of free radicals, damaging to tissues.

The T allele affects proteolytic cleavage of the enzyme thereby reducing nitric oxide bio-availability in the blood vessel wall and promoting atherosclerosis, as a result it is associated with atherosclerosis, essential hypertension, end-stage renal disease and pre-eclampsia.⁵⁹⁻⁶¹ The T allele has also been associated with a number of cancers.^{62,63}

GT and TT genotype individuals might show greater beneficial effects of n-3 PUFA consumption.⁶⁰

3. Genetic variants associated with athletic performance

It is now two decades since the first gene was identified as an indicator of sporting performance.⁶⁴ This gene, ACE, has now been extensively studied and been shown to be predictive of endurance and power-oriented performance. Numerous other genes have also been studied and fourteen genetic polymorphisms have been included in the DNA Sport test. They have been identified as those with high scientific credibility and those which can make a difference to the training and lifestyle intervention of an athlete, based on knowledge of his/her genes.

3.1. Angiotensinogen (AGT): Met238Thr (T>C)

The angiotensinogen (AGT) protein, produced mainly by the liver, is an important component of the renin-angiotensin system. AGT is cleaved by renin to form angiotensin I, the precursor of angiotensin II that regulates vascular resistance and sodium homeostasis through vasoconstriction, determining blood pressure. High levels of AGT can lead to an increase in the production of angiotensin II resulting in hypertension.⁶⁵

The Met238Thr (T>C) polymorphism within AGT results in a threonine amino acid, instead of a methionine, at position 238 of the protein. This variant has been associated with higher plasma levels of AGT, blood pressure regulation and the risk of hypertension-associated disorders.²⁷ AGT Met238Thr (T>C) has also been associated with elite power and strength performance.

When comparing the frequency of the AGT Met238Thr (T>C) variation between top endurance and power Spanish athletes as well as controls, Gomez-Gallego *et al.* (2009) showed that the CC genotype was significantly over-represented in power athletes. This indicated that the CC genotype may favour power sports performance, possibly due to the higher activity of angiotensin II that acts as a skeletal muscle growth factor.⁶⁶ Zarębska *et al.* (2013) repeated this study in Polish athletes and found that the frequency of the CC genotype was significantly higher in power athletes in comparison to control and endurance groups.⁶⁷ This AGT polymorphism has recently been investigated with regards to response to a 12 week aerobics program.⁶⁸ Participants completed physical power tests before and after the intervention, such as sprint tests and jump tests. All power-related variables improved significantly in response to the aerobic training, however significantly greater improvements were seen in the jump-based power variables in C allele carriers in comparison to TT homozygotes.⁶⁸ This indicates that the AGT gene may influence power and strength improvements in response to training.

The C allele has been associated with significantly higher levels of plasma AGT and an increased resting blood pressure.⁶⁹ The Heritage Family Study investigated the AGT Met238Thr polymorphism and blood pressure response to an endurance exercise training program. T allele carriers showed significantly greater decreases in diastolic blood pressure in comparison to CC homozygotes.⁷⁰ Additionally, a longitudinal study demonstrated that regular moderate-intensity exercise attenuates age-related increases in systolic blood pressure in TT homozygotes.⁷¹

AGT CC homozygotes are more susceptible to left ventricular cardiac hypertrophy with long term exercise training,⁷² and even in response to a 17 week training program.⁷³ The AGT CC and ACE DD genotypes have also been shown to interact in the effect of exercise training on cardiac mass.⁷⁴

3.2. Angiotensin I-Converting Enzyme (ACE): Insertion (I)/Deletion (D)

The angiotensin I converting enzyme (ACE) functions in regulating circulatory homeostasis, through the synthesis of the vasoconstrictor angiotensin II. ACE hydrolyses angiotensin I to angiotensin II, which also drives aldosterone synthesis, and the degradation of bradykinin, a potent vasodilator.

A polymorphism exists in the ACE gene that results in the insertion (I allele) or the absence (D allele) of a 287 base pair fragment. Inter-individual variation in ACE plasma levels can vary as much as a 5-fold, approximately half of this variation is thought to be due to the ACE I/D polymorphism.⁷⁵

The D allele is associated with increased activity of the ACE enzyme, likely causing increased production of the angiotensin II vasoconstrictor and degradation of the vasodilator bradykinin. Higher circulating ACE activity has been significantly associated with baseline quadriceps muscle strength.⁷⁶ The I allele on the other hand has been associated with lower ACE activity.

The ACE DD genotype has been associated with greater lean body mass, muscle volume, and percentage of Type II fast twitch muscle fibres.^{77,78} Various studies have shown that the D allele is over-represented in power and strength athletes.

The D allele or DD genotype has shown to be over-represented in short and middle distance swimmers.⁷⁹⁻⁸¹ Nazarov *et al.* (2001) showed a significantly higher frequency of the ACE D allele in short distance athletes (<1 minute),⁸² similarly shown by Papadimitriou *et al.* (2009) in track and field sprinters.⁸³ Recently, the ACE D allele has been significantly associated with greater standing log jump results.⁸⁴ Giaccaglia *et al.* (2008) found following an 18-month walking and light weight training exercise programme involving elderly, obese individuals, those with the DD genotype had greater gains in knee extensor muscle strength in comparison to subjects with the II genotypes.⁸⁵ Following a 12-week high speed training program in older woman, Pereira *et al.* (2013) showed that D allele carriers had greater improvements in maximum strength and power.⁸⁶

Individuals with the ACE DD genotype should be aware of an increased risk left ventricular hypertrophy as a result of exercise training.⁶⁴

The I allele on the other hand has been linked to aerobic endurance athletic ability and a greater percentage of type I slow twitch muscle fibres,⁷⁸ higher VO₂max, greater aerobic work efficiency, improved fatigue resistance, higher peripheral tissue oxygenation during exercise, greater aerobic power response to training, greater cardiac output and maximal power output in athletes.⁸⁷

Myerson *et al.* (1999) showed that the frequency of the I allele increased with increasing running distance in Olympic runners, revealing that the I allele may be associated with improved endurance performance.⁷⁹ The ACE I allele has also been linked to elite endurance rowing,^{88,89} and shown to be

over-represented in elite cyclists, long distance runners and handball players.⁹⁰ Collins *et al.* (2004) noted an association between the I allele and endurance performance of fastest triathlon finishers.⁹¹ Similarly, Shenoy *et al.* (2010) showed a greater frequency of I alleles in triathletes compared to controls.⁹² The II genotype is also over-represented in endurance runners,^{93,94} with frequency of the allele increasing with running distance.⁹⁵ Znazen *et al.* (2016) recently noted that the ACE I/D polymorphism frequency was significantly different between power and endurance athletes, suggesting that this may play an important role in sporting adaptation to physical training.⁹⁶ A recent meta-analysis by Ma *et al.* (2013) further solidified the role of the II genotype in endurance athlete status.⁹⁷

A 2016 multi-cohort study examined the ACE I/D polymorphism and sprint time in elite athletes. ACE DD genotype carriers had faster 200m and 400m sprint times compared to those with the II genotype. The study found that the ACE D allele, in a dominant model, accounted 1.48 % of sprint time variance. Although multiple genetic variants and environmental factors influence sprint time, a combination of the ACE and ACTN3 gene variants may be the difference between a world record and only making the Olympic final.⁹⁸

3.3. Bradykinin Receptor B2 (BDKRB2): Exon 1 +9/-9 (C>T)

Bradykinin is an endothelium-dependent vasodilator and acts via the Bradykinin B2 Receptor (BDKRB2). ACE is able to degrade the vasodilator bradykinin by forming a complex with BDKRB2, this is thought to play a role in the cross-talk between the renin-angiotensin and the kallikrein-kinin systems.

The absence (-9) of a 9 base pair repeat in exon 1 of BDKRB2 is associated with higher gene transcriptional activity and mRNA expression of BDKRB2 (T allele).⁹⁹ Levels of bradykinin may effect skeletal muscle glucose uptake and blood flow during exercise. The BDKRB2 -9 (T) allele has been positively associated with endurance exercise performance, likely due to enhanced skeletal muscle mechanical and metabolic efficiency (the amount of energy used per unit of power output) with this allele.⁹⁹

In 2004, Williams *et al.* showed that the BDKRB2 -9 (T) allele is associated with greater muscle contraction efficiency in healthy individuals as well as Olympic track athletes.⁹⁹ This study also showed that a combination of the ACE I allele and BDKRB2 -9 (T) allele yields an additive effect, significantly associated with aerobic endurance capacity.

Saunders *et al.* 2006, showed that the BDKRB2 -9-9 (TT) genotype occurred at a significantly higher frequency in an ironman triathlete group compared to controls. The -9 (T) allele also appears to be associated with a faster finishing time.¹⁰⁰ A higher frequency of the BDKRB2 TT genotype has also been noted in marathon runners.¹⁰¹

The BDKRB2 +9+9 (CC) genotype has been associated with greater thirst and fluid weight loss during an Ironman triathlon event; this would not be advantageous for endurance activities.¹⁰²

A subjects' BDKRB2 +9/-9 (C>T) genotype should be interpreted with caution as some studies have failed to replicate the association of this polymorphism with endurance athlete status.^{103,104}

3.4. Vascular Endothelial Growth Factor (VEGF): -634 C>G

Vascular Endothelial Growth Factor (VEGF) is an endothelial cell proliferator involved in blood vessel growth. Angiogenesis is a key factor in the physiological adaptation to aerobic exercise, required to match the needs of the active tissue. Levels of VEGF affect peripheral circulation and capillary blood flow, muscle oxygenation and VO₂max (maximal oxygen consumption).¹⁰⁵ VEGF promotes angiogenesis in response to aerobic exercise, an important contributor to training-induced changes in VO₂max and increases in capillary-muscle fibre ratio and capillary density.¹⁰⁶ Coffey *et al.* (2006) noted an approximate 4-fold increase in muscle VEGF levels in response to both resistance and endurance exercise training sessions,¹⁰⁷ while Gustafsson *et al.* (2005) showed a 178% increase in muscle VEGF mRNA expression following dynamic resistance exercise.¹⁰⁸

The -634 C>G polymorphism in the promoter region of the VEGF gene has been linked to altered expression of the VEGF protein.¹⁰⁹ This gene variant has been associated with aerobic capacity and endurance athlete status.¹⁰⁵

The -634 C>G C allele is associated with greater expression of VEGF and a higher VO₂max level at baseline and in response to aerobic training.¹⁰⁵ Ahmetov *et al.* (2008) noted that the frequency of the VEGF C allele is significantly higher in endurance athletes in comparison to control subjects.¹¹⁰ This was replicated in another large cohort of 1423 athletes and 1132 controls.¹¹¹ A correlation is also seen between the VEGF C allele and higher levels of aerobic performance (based on maximum power and VO₂max).¹¹⁰ The combination of studies suggest that increased VEGF expression, as seen with the C allele, is associated with greater muscle efficiency, advantageous for aerobic performance.

3.5. Nuclear Respiratory Factor 2 (NRF-2); GABPB1: A>G

The GA binding protein transcription factor β subunit 1 (GABPB1), also known as nuclear respiratory factor 2 (NRF2) is one of the two nuclear respiratory factors, which induce mitochondrial biogenesis. The NRF2 transcription factor regulates genes involved in mitochondrial function, influencing respiratory capacity and the rate of ATP production during exercise. The stimulation of mitochondrial biogenesis by exercise is dependent on increased NRF2 gene expression. NRF2 also regulates mitochondrial transcription factor A (TFAM), cytochrome c, and heme biosynthesis proteins.¹¹²

An A>G variant within NRF2 (GABPB1) has been linked to VO₂max and endurance capacity. In 2007, He *et al.* showed that a greater training response to an 18-wk endurance training program was seen in NRF2 rs7181866 G allele carriers. When three SNPs on the NRF-2 gene (rs12594956, rs8031031 and rs7181866) those with the ATG haplotype had a 57% higher training response, determined by running

economy and VO₂max, compared non-carriers of the particular alleles.¹¹² Eynon *et al.* (2009) found a significantly higher proportion of the AG genotype in endurance athletes, compared to sprinters, and control subjects. This association was even more pronounced when looking at elite international athletes, suggesting an association between the G allele and endurance athlete status.¹¹³ Maciejewska-Karłowska *et al.* (2012) observed that the frequency of the AG genotype was significantly higher in a group of rowers compared to controls.¹¹⁴

The collective studies support the notion that genetic variation within the NRF2 gene belongs to group of polymorphisms that are associated with endurance performance.

3.6. Peroxisome Proliferator-Activated Receptor-Gamma Coactivator-1 (PPARGC1A): Gly482Ser (G>A)

Peroxisome Proliferator-Activated Receptor-Gamma Coactivator-1 (PPARGC1A), otherwise known as PGC-1 α , is a co-activator of the PPAR family of proteins and plays an essential role in mitochondrial function and biogenesis, fatty acid oxidation, energy homeostasis and glucose utilisation, thermogenesis, adipocyte differentiation, angiogenesis, and Type I slow twitch muscle fibre specialization.¹¹⁵ PPARGC1A is expressed in tissues with high energy demands, where a high mitochondrial density is required. Exercise training results in increased expression of PPARGC1A (preferentially in Type 1 muscle fibres), a key factor in exercise-induced mitochondrial biogenesis.¹¹⁶

A Gly482Ser (G>A) polymorphism exists in PPARGC1A that influences expression of this gene; the A allele resulting in decreased levels of PPARGC1A mRNA,¹¹⁵ consequently associated with reduced aerobic improvements from exercise training. The GG genotype on the other hand has been associated with a greater VO₂max, reflective of aerobic capacity.¹¹⁵ Interestingly, the frequency of the A allele has been observed to be higher in individuals with Type 2 diabetes.

Lucia *et al.* (2005) showed that the frequency of the unfavourable A allele was significantly lower in world class Spanish endurance athletes in comparison to controls.¹¹⁷ Individuals with the G allele also had greater improvements in aerobic fitness following a 9 month intervention study in comparison to A allele carriers.¹¹⁸ Akhmetov *et al.* (2007) observed a higher frequency of the PPARGC1A G allele in athletes in comparison to controls; the G allele also being associated high values of aerobic performance.¹¹⁹ The high frequency of the G allele in endurance athletes was seen again in a large study combining multiple genetic variants linked to endurance athlete status.¹¹¹ This has similarly been repeated in Israeli,¹²⁰ Polish and Russian¹²¹ endurance athletes.

Conversely, Gineviciene *et al.* (2016) recently found an association between the AA genotype and strength/power athlete status, with the AA genotype being particularly favourable for powerlifters.¹²²

3.7. Peroxisome Proliferator-Activated Receptor-Alpha (PPARA): Intron 7 G>C

Peroxisome Proliferator-Activated Receptor-Alpha (PPARA) transcription factor is a regulator of lipid and glucose metabolism and energy homeostasis in the liver, skeletal muscle and heart. Activation of PPARA promotes the uptake, utilisation, and oxidation of fatty acids and plays an important role in the adaptive response to endurance training. PPARA expression levels are higher in type I (slow-twitch) muscle fibres than in type II (fast-twitch) muscle fibres.¹²³ Endurance training increases the use of non-plasma fatty acids in energy homeostasis and may enhance skeletal muscle oxidative capacity by influencing PPARA gene expression.

A G>C variation exists in intron 7 of the PPARA gene. The G allele leads to increased expression of PPARA and has been associated with endurance performance. Conversely the C allele is linked to anaerobic (power) performance predominantly reliant on muscle glucose metabolism.¹²³

Ahmetov *et al.* (2006) noted a higher frequency of the PPARA GG genotype in endurance athletes, compared to an increasingly higher frequency of the C allele with increasing anaerobic component of physical performance. Individuals with the GG genotype also had a significantly higher percentage of slow twitch muscle fibres in comparison to CC homozygotes.¹²³ The same authors found that the frequency of the PPARA G allele was significantly greater in elite rowers, and associated with higher values of aerobic performance, than in control subjects.¹¹⁹ Maciejewska *et al.* (2011) similarly noted a higher frequency of the G in elite rowers compared to sedentary controls.

The PPARA C allele has been associated with exercise-induced left ventricular growth. Jamshidi *et al.* (2002) noted that individuals with the C allele had a significantly greater increase in left ventricular mass following a 10-week training program in comparison to GG homozygotes.¹²⁴

3.8. Beta 2 Adrenergic Receptor (ADRB2): Arg16Gly (A>G); Gln27Glu (C>G)

The adrenergic receptors are a family regulatory proteins that mediate the physiological effects of the epinephrine and the norepinephrine. ADRB2 encodes the β_2 adrenergic receptor which plays a part in the regulation of the cardiac, pulmonary, vascular, endocrine and the central nervous system.¹²⁵ Adrenaline acts predominantly via ADRB2 in order to maintain blood glucose levels during prolonged exercise, by promoting glycogenolysis.

A variant that causes an arginine to glycine change at position 16 of ADRB2 (A>G) has been associated with altered vasodilator responses to catecholamines during stress, and so modulates the pressor response (increasing cardiac output) to exercise, as driven by the sympathetic nervous system. A glutamine to glutamate polymorphism at position 27 of the gene (C>G) has been associated with endurance performance and the ability to lose weight with exercise training.

Kochanska-Dziurawicz *et al.* (2013) measured plasma concentrations of adrenaline, noradrenaline and ADRB2 gene expression before and after exercise in trained ice hockey players. Performed work by the athletes was found to be dependent on initial levels of noradrenaline in plasma and ADRB2 mRNA in peripheral blood mononuclear cells.¹²⁵ It has been suggested that the endogenous catecholamine stimulation of ADRB2 that occurs with heavy exercise may result in enhanced cardiovascular function in individuals with one particular genotype over another.¹²⁶

Arg16Gly (A>G)

Snyder *et al.* (2006) found that individuals with Arg (A) allele had lower receptor density and resting cardiac output,¹²⁶ and suggest a correlation between these factors. Eisenach *et al.* (2004) noted a greater heart rate response to exercise in individuals with the Arg/Arg (AA) genotype, which may aid in maintaining an increased cardiac output during peripheral vasodilation (decreased total peripheral resistance).¹²⁷ In an investigation of the Arg16Gly variant in athletes, Wolfarth *et al.* showed that the 16Arg (A) allele was significantly overrepresented in elite endurance athletes compared to sedentary controls, suggesting a positive association between the Arg16Gly polymorphism and endurance performance status.¹²⁸ Tsianos *et al.* (2010) reported that marathon runners with the ADRB2 Arg (A) allele had significantly faster running times compared to (Gly) G allele carriers, implicating this variant in endurance running. Performance.¹⁰¹ The Arg (A) has also been associated with greater peak VO₂max in heart failure patients.¹²⁹

Gln27Glu (C>G)

Moore *et al.* (2001) investigated the ADRB2 Gln27Glu variant in female elite endurance athletes post menopause. Individuals with the Gln (C) allele had a greater VO₂max than Glu/Glu (GG) homozygotes and were subsequently associated with elite endurance performance.¹³⁰ When looking at the variants in the ADRB2 gene in athletes of varying sporting demands (both power/strength and endurance), Sawczuk *et al.* (2013) found that the Arg 16Gly G allele and the Gln7Glu G allele, were overrepresented in strength/power athletes compared to controls.¹³¹

The Glu/Glu (GG) genotype has also been associated with greater fat deposition and obesity. Macho-Azcarate *et al.* 2003 noted that women with the Glu/Glu (GG) genotype had reduced fat oxidation in response to exercise in comparison to those with the Gln (C) allele.¹³²

3.9. Thyrotropin Releasing Hormone Receptor (TRHR): rs7832552 C>T

The TRHR gene encodes the Thyrotropin-releasing hormone (TRH) receptor. TRH stimulates the release of thyroid-stimulating hormone (TSH) and prolactin from the anterior pituitary. TSH in turn stimulates the release of Thyroxin (T4) and Triiodothyronine (T3) from the thyroid gland. Thyroxin is necessary for the anabolic action of the growth factor GH-IGF1, essential in muscle protein balance and adaptive changes to load.¹³³ T3 and T4 are also involved in increasing metabolic rate and aiding catecholamines in the mobilization of fuels during exercise.

Two gene variants within TRHR (rs7832552 C>T and rs16892496 T>G), have been significantly associated with lean body mass (skeletal muscle mass).¹³³ These variants are in strong linkage disequilibrium and the minor T and G alleles have respectively been found to be favourable for increased lean body mass and therefore strength and power activities.

Liu *et al.* (2009) conducted a genome wide association study (GWAS) on lean body mass variation. TRHR was identified as an important genetic determinant of lean body mass. Subjects with the TT rs7832552 genotype had 2.55 kg greater lean body mass compared to those with the C allele.¹³³ CC This same association was observed in three replication studies, totalling more than 6000 participants, and meta analyses of the initial GWAS study and replication studies achieved significant association.

3.10. Alpha-actinin 3 (ACTN3): 577 R>X

Alpha-actinins are a family of actin-binding proteins that maintain the cytoskeleton. The ACTN3 gene encodes alpha-actinin 3, which is only present in type II (fast) muscle fibres and also has a low expression in brain tissue. Alpha-actinins are found at the Z-line of the muscle where they anchor actin filaments and help to maintain the structure of the sarcomere. Alpha-actinins also interact with a signalling factor that plays a role in the specialisation of muscle fibre type, diameter and metabolism.

The ACTN3 gene contains a polymorphism which results in two versions of ACTN3: a functional R allele and a null X allele. The genotype that is homozygous for the X allele (XX) is completely deficient in alpha-actinin 3. The percentage of people with two X alleles range from about 1% in African populations to about 18% in Europeans and 25% in East Asians.

It has been demonstrated that male and female elite Sprint athletes have significantly higher frequencies of the R alleles than endurance athletes and sedentary controls.^{134,135} Additionally, when Roth *et al* (2008) studied elite and local black and white bodybuilders, they showed a significantly lower XX frequency in the bodybuilders compared to controls.¹³⁶ None of the black bodybuilders were homozygous for the X allele. In a cycling study, it was found that subjects who were not ACTN3 deficient had greater peak power outputs and ventilatory thresholds than cyclists with the XX genotype.¹³⁷

When studying the ACTN3 gene for endurance performance, mixed results have been found. Shang *et al* (2010) showed that the XX genotype was significantly over-represented in the female endurance athletes compared to controls.¹³⁸ However, no relationship between genotype and endurance-status was found in the male athletes. Doring *et al* (2010) found that the XX genotype was similar in prevalence between 316 endurance athletes and 304 controls.¹³⁹ Additionally, Saunders *et al* (2007) found when studying 457 male Ironman triathletes compared with 243 controls, that the XX polymorphism was not associated with ultra-endurance performance.¹⁴⁰ These observations have also been replicated in professional cyclists and Olympic endurance runners.¹⁴¹ Collectively these studies suggest that the XX genotype has not yet been shown to be a good marker for endurance performance.

3.11. Vitamin D Receptor (VDR): Taq1 T>C

The Vitamin D Receptor (VDR) gene encodes the nuclear hormone receptor for vitamin D3 and is expressed by cells in most organs. VDR activation is essential for the maintenance of calcium and phosphorus levels in the blood and bone. VDR regulates the expression of many genes involved in essential bodily functions and has been identified as a regulator of skeletal muscle;¹⁴² VDR binds vitamin D metabolites in muscle cells.

Approximately 200 variants exist in the VDR gene. The VDR Taq1 polymorphism is a T>C (T>t) change; this variant has been shown to be in linkage disequilibrium with the VDR Bsm1 G>A (b>B) variant. VDR Taq1 is reported in DNA Sport, however VDR Bsm1 is also commonly referred to in the literature. The VDR Taq1 CC (commonly called tt) and VDR Bsm1 AA (commonly called BB) homozygotes have been associated with greater muscle strength.

Geusens *et al.* (1997) investigated VDR polymorphisms and muscle strength. Quadriceps strength was found to be significantly higher in non-obese elderly women with the VDR Bsm1 bb genotype (Taq1 TT).¹⁴³ In contrast Grundberg *et al.* (2004) found that healthy women with the VDR Bsm1 BB (Taq1 CC) genotype, had greater hamstring isokinetic strength than the bb (TT) genotype group.¹⁴⁴ Additionally, Wang *et al.* (2006) showed that VDR Bsm1 bb (Taq1 TT) genotypes demonstrated lower knee flexion concentric torque than the B allele (Taq1 C allele) carriers.¹⁴⁵ Windelinckx *et al.* (2007) demonstrated that men with the Bsm1 AA genotype (BB) and Taq1 CC genotype (tt) had higher isometric and concentric quadriceps strength in comparison to their counterparts.¹⁴⁶ Bahat *et al.* (2010) similarly replicated the association of the VDR Bsm1 BB (Taq1 CC) genotype with greater extensor strength of the knee.¹⁴⁷

Vitamin D supplementation may improve muscle strength through an increase in the size and amount of type II (fast twitch) muscle fibres in individuals with low vitamin D status.¹⁴⁸ Sufficient Vitamin D status is important in athletes and has also been linked to the prevention of stress fractures.¹⁴² The VDR Bsm1 BB genotype (Taq1 CC) has additionally been associated with lower bone mineral density.¹⁴⁹

3.12. Cytochrome P450 1A2 (CYP1A2): -163 C>A

Caffeine is a central nervous system and metabolic stimulant that is used to reduce physical fatigue. In athletics, moderate doses of caffeine have been known to improve both sprint and endurance performance.¹⁵⁰ Cytochrome P450 1A2 (CYP1A2) is one of the main enzymes that catalyse the oxidation of caffeine in humans.

A C>A variant exists in intron 1 of CYP1A2, affecting caffeine metabolism. Individuals with the C allele are associated with a reduced ability to metabolise caffeine.¹⁵¹

Womack *et al.* (2012) investigated the ergogenic effects of caffeine and the CYP1A2 polymorphism in trained cyclists. Individuals with the AA genotype experienced a greater ergogenic effect and

enhanced performance benefit with caffeine ingestion, in comparison to C allele carriers. This study was performed using 6mg/kg of caffeine 1 hour prior to performance.¹⁵⁰ Pataky *et al.* (2016) further demonstrated that cyclists with the CYP1A2 A allele benefited from caffeine ingestion with or without a caffeine mouth rinse 1 hour prior to exercise.¹⁵²

Conversely, Salinero *et al.* (2017) found no significant effect of the CYP1A2 genotype on performance with 3mg/kg caffeine intake 1 hour prior to exercise.¹⁵³

Individuals with the CYP1A2 AA genotype are likely to metabolise and clear caffeine at a faster rate than those with the C allele. C allele carriers should experiment with taking caffeine in more than 1 hour before the start of a race or event in order to gain from the ergogenic effects. Individuals with the CYP1A2 AA genotype are likely to benefit from caffeine intake \leq 1 hour before a race or event and depending on the length of the event, consider an additional dose.

3.13. Circadian Locomotor Output Cycles Kaput (CLOCK): 3111 T>C

Multiple variables influence sports or exercise performance; more recently, circadian rhythms and chronotype have been noted as considerable factors. Chronotype refers to morning preference, evening preference or neither, and affects several physiological variables. An individual's chronotype may affect the time of day they are most likely to be at peak performance.^{154,155}

A systematic review by Vitale *et al.* (2017) on chronotype and exercise concluded that morning and evening preference effects athletic performance.¹⁵⁵ Individuals that have a morning preference chronotype perceived less effort when performing a submaximal physical task in the morning in comparison to individuals with evening preference or neither. Additionally, individuals with morning preference achieved better race times in morning events in comparison to those with evening preference or neither.

Circadian Locomotor Output Cycles Kaput (CLOCK) is an essential element of the persistence and period of circadian rhythms, and effects chronotype and sleep. The C allele of the CLOCK 311 T>C variant has been associated with evening preference and delayed sleep timing.¹⁵⁶⁻¹⁵⁸ This gene has also been linked to weight loss and eating behaviour.¹⁵⁷

Genetic variants that underlie chronotype (morning or evening preference) and the effect on performance, need to be further studied in order to confirm previous findings. Coaches and trainers should however consider the chronotype effect, habitual time of day of training sessions, and time of day of competition in the lead up to an event.¹⁵⁹

4. References

1. Kannus, P. Structure of the tendon connective tissue. *Scand. J. Med. Sci. Sports* **10**, 312–320 (2000).
2. Khoschnau, S. *et al.* Type I collagen alpha1 Sp1 polymorphism and the risk of cruciate ligament ruptures or shoulder dislocations. *Am. J. Sports Med.* **36**, 2432–2436 (2008).
3. Posthumus, M. *et al.* Genetic risk factors for anterior cruciate ligament ruptures: COL1A1 gene variant. *Br. J. Sports Med.* **43**, 352–356 (2009).
4. Posthumus, M., September, A. V., Schwellnus, M. P. & Collins, M. Investigation of the Sp1-binding site polymorphism within the COL1A1 gene in participants with Achilles tendon injuries and controls. *J. Sci. Med. Sport Sports Med. Aust.* **12**, 184–189 (2009).
5. Collins, M., Posthumus, M. & Schwellnus, M. P. The COL1A1 gene and acute soft tissue ruptures. *Br. J. Sports Med.* **44**, 1063–1064 (2010).
6. Mann, V. *et al.* A COL1A1 Sp1 binding site polymorphism predisposes to osteoporotic fracture by affecting bone density and quality. *J. Clin. Invest.* **107**, 899–907 (2001).
7. Wenstrup, R. J. *et al.* Type V Collagen Controls the Initiation of Collagen Fibril Assembly. *J. Biol. Chem.* **279**, 53331–53337 (2004).
8. Collins, M. & Posthumus, M. Type V collagen genotype and exercise-related phenotype relationships: a novel hypothesis. *Exerc. Sport Sci. Rev.* **39**, 191–198 (2011).
9. Mokone, G. G., Schwellnus, M. P., Noakes, T. D. & Collins, M. The COL5A1 gene and Achilles tendon pathology. *Scand. J. Med. Sci. Sports* **16**, 19–26 (2006).
10. September, A. V. *et al.* Variants within the COL5A1 gene are associated with Achilles tendinopathy in two populations. *Br. J. Sports Med.* **43**, 357–365 (2009).
11. Posthumus, M. *et al.* The COL5A1 gene is associated with increased risk of anterior cruciate ligament ruptures in female participants. *Am. J. Sports Med.* **37**, 2234–2240 (2009).
12. Altinisik, J. *et al.* The BstUI and DpnII Variants of the COL5A1 Gene Are Associated With Tennis Elbow. *Am. J. Sports Med.* **43**, 1784–1789 (2015).
13. Brown, K. L. *et al.* Polymorphisms within the COL5A1 gene and regulators of the extracellular matrix modify the risk of Achilles tendon pathology in a British case-control study. *J. Sports Sci.* **35**, 1475–1483 (2017).
14. Laguette, M.-J., Abrahams, Y., Prince, S. & Collins, M. Sequence variants within the 3'-UTR of the COL5A1 gene alters mRNA stability: implications for musculoskeletal soft tissue injuries. *Matrix Biol. J. Int. Soc. Matrix Biol.* **30**, 338–345 (2011).

15. Dressler, M. R. *et al.* A potential mechanism for age-related declines in patellar tendon biomechanics. *J. Orthop. Res. Off. Publ. Orthop. Res. Soc.* **20**, 1315–1322 (2002).
16. Collins, M., Mokone, G. G., September, A. V., van der Merwe, L. & SchwelInus, M. P. The COL5A1 genotype is associated with range of motion measurements. *Scand. J. Med. Sci. Sports* **19**, 803–810 (2009).
17. Brown, J. C., Miller, C.-J., SchwelInus, M. P. & Collins, M. Range of motion measurements diverge with increasing age for COL5A1 genotypes. *Scand. J. Med. Sci. Sports* **21**, e266-272 (2011).
18. Posthumus, M., SchwelInus, M. P. & Collins, M. The COL5A1 gene: a novel marker of endurance running performance. *Med. Sci. Sports Exerc.* **43**, 584–589 (2011).
19. Brown, J. C., Miller, C.-J., Posthumus, M., SchwelInus, M. P. & Collins, M. The COL5A1 gene, ultra-marathon running performance, and range of motion. *Int. J. Sports Physiol. Perform.* **6**, 485–496 (2011).
20. Mikic, B., Schalet, B. J., Clark, R. T., Gaschen, V. & Hunziker, E. B. GDF-5 deficiency in mice alters the ultrastructure, mechanical properties and composition of the Achilles tendon. *J. Orthop. Res. Off. Publ. Orthop. Res. Soc.* **19**, 365–371 (2001).
21. Tashiro, T. *et al.* Effect of GDF-5 on ligament healing. *J. Orthop. Res. Off. Publ. Orthop. Res. Soc.* **24**, 71–79 (2006).
22. Ozasa, Y. *et al.* A comparative study of the effects of growth and differentiation factor 5 on muscle-derived stem cells and bone marrow stromal cells in an in vitro tendon healing model. *J. Hand Surg.* **39**, 1706–1713 (2014).
23. Posthumus, M. *et al.* Components of the transforming growth factor-beta family and the pathogenesis of human Achilles tendon pathology--a genetic association study. *Rheumatol. Oxf. Engl.* **49**, 2090–2097 (2010).
24. Egli, R. J. *et al.* Functional analysis of the osteoarthritis susceptibility-associated GDF5 regulatory polymorphism. *Arthritis Rheum.* **60**, 2055–2064 (2009).
25. Vaes, R. B. A. *et al.* Genetic variation in the GDF5 region is associated with osteoarthritis, height, hip axis length and fracture risk: the Rotterdam study. *Ann. Rheum. Dis.* **68**, 1754–1760 (2009).
26. Miyamoto, Y. *et al.* A functional polymorphism in the 5' UTR of GDF5 is associated with susceptibility to osteoarthritis. *Nat. Genet.* **39**, 529–533 (2007).
27. Southam, L. *et al.* An SNP in the 5'-UTR of GDF5 is associated with osteoarthritis susceptibility in Europeans and with in vivo differences in allelic expression in articular cartilage. *Hum. Mol. Genet.* **16**, 2226–2232 (2007).
28. Huuskonen, A. *et al.* A Common Variation in the Promoter Region of Interleukin-6 Gene Shows Association with Exercise Performance. *J. Sports Sci. Med.* **8**, 271–277 (2009).
29. Yamin, C. *et al.* IL6 (-174) and TNFA (-308) promoter polymorphisms are associated with systemic creatine kinase response to eccentric exercise. *Eur. J. Appl. Physiol.* **104**, 579–586 (2008).
30. Buxens, A. *et al.* Can we predict top-level sports performance in power vs endurance events? A genetic approach. *Scand. J. Med. Sci. Sports* **21**, 570–579 (2011).

31. Ruiz, J. R. *et al.* Can we identify a power-oriented polygenic profile? *J. Appl. Physiol.* **108**, 561–566 (2010).
32. Robson, P. Elucidating the unexplained underperformance syndrome in endurance athletes: the interleukin-6 hypothesis. *Sports Med. Auckl. NZ* **33**, 771–781 (2003).
33. Robson-Ansley, P. *et al.* The effect of repeated endurance exercise on IL-6 and sIL-6R and their relationship with sensations of fatigue at rest. *Cytokine* **45**, 111–116 (2009).
34. McKenzie, J. A. *et al.* Influence of the interleukin-6 -174 G/C gene polymorphism on exercise training-induced changes in glucose tolerance indexes. *J. Appl. Physiol. Bethesda Md 1985* **97**, 1338–1342 (2004).
35. Rose-John, S., Waetzig, G. H., Scheller, J., Grötzing, J. & Seeger, D. The IL-6/sIL-6R complex as a novel target for therapeutic approaches. *Expert Opin. Ther. Targets* **11**, 613–624 (2007).
36. Gray, S. R. *et al.* The response of circulating levels of the interleukin-6/interleukin-6 receptor complex to exercise in young men. *Cytokine* **47**, 98–102 (2009).
37. Galicia, J. C. *et al.* Polymorphisms in the IL-6 receptor (IL-6R) gene: strong evidence that serum levels of soluble IL-6R are genetically influenced. *Genes Immun.* **5**, 513–516 (2004).
38. Reich, D. *et al.* Admixture mapping of an allele affecting interleukin 6 soluble receptor and interleukin 6 levels. *Am. J. Hum. Genet.* **80**, 716–726 (2007).
39. Gray, S. R., Robinson, M. & Nimmo, M. A. Response of plasma IL-6 and its soluble receptors during submaximal exercise to fatigue in sedentary middle-aged men. *Cell Stress Chaperones* **13**, 247–251 (2008).
40. Obisesan, T. O. *et al.* C-reactive protein genotypes affect baseline, but not exercise training-induced changes, in C-reactive protein levels. *Arterioscler. Thromb. Vasc. Biol.* **24**, 1874–1879 (2004).
41. Eiriksdottir, G. *et al.* The interaction of adiposity with the CRP gene affects CRP levels: age, gene/environment susceptibility-Reykjavik study. *Int. J. Obes.* **2005** **33**, 267–272 (2009).
42. Kullo, I. J., Khaleghi, M. & Hensrud, D. D. Markers of inflammation are inversely associated with VO₂ max in asymptomatic men. *J. Appl. Physiol. Bethesda Md 1985* **102**, 1374–1379 (2007).
43. Kuo, H.-K., Yen, C.-J., Chen, J.-H., Yu, Y.-H. & Bean, J. F. Association of cardiorespiratory fitness and levels of C-reactive protein: data from the National Health and Nutrition Examination Survey 1999-2002. *Int. J. Cardiol.* **114**, 28–33 (2007).
44. Neubauer, O., König, D. & Wagner, K.-H. Recovery after an Ironman triathlon: sustained inflammatory responses and muscular stress. *Eur. J. Appl. Physiol.* **104**, 417–426 (2008).
45. Depner, C. M., Kirwan, R. D., Frederickson, S. J. & Miles, M. P. Enhanced inflammation with high carbohydrate intake during recovery from eccentric exercise. *Eur. J. Appl. Physiol.* **109**, 1067–1076 (2010).
46. Ross, M. L. R. *et al.* Cytokine responses to carbohydrate ingestion during recovery from exercise-induced muscle injury. *J. Interferon Cytokine Res. Off. J. Int. Soc. Interferon Cytokine Res.* **30**, 329–337 (2010).

47. Robson-Ansley, P., Walshe, I. & Ward, D. The effect of carbohydrate ingestion on plasma interleukin-6, hepcidin and iron concentrations following prolonged exercise. *Cytokine* **53**, 196–200 (2011).
48. Phillips, T., Childs, A. C., Dreon, D. M., Phinney, S. & Leeuwenburgh, C. A dietary supplement attenuates IL-6 and CRP after eccentric exercise in untrained males. *Med. Sci. Sports Exerc.* **35**, 2032–2037 (2003).
49. McGinley, C., Shafat, A. & Donnelly, A. E. Does antioxidant vitamin supplementation protect against muscle damage? *Sports Med. Auckl. NZ* **39**, 1011–1032 (2009).
50. Moldoveanu, A. I., Shephard, R. J. & Shek, P. N. Exercise elevates plasma levels but not gene expression of IL-1beta, IL-6, and TNF-alpha in blood mononuclear cells. *J. Appl. Physiol. Bethesda Md* **1985** **89**, 1499–1504 (2000).
51. Lakka, H.-M. *et al.* The TNF-alpha G-308A polymorphism is associated with C-reactive protein levels: the HERITAGE Family Study. *Vascul. Pharmacol.* **44**, 377–383 (2006).
52. Nicklas, B. J. *et al.* Physical function and its response to exercise: associations with cytokine gene variation in older adults with knee osteoarthritis. *J. Gerontol. A. Biol. Sci. Med. Sci.* **60**, 1292–1298 (2005).
53. Nieman, D. C. *et al.* Muscle cytokine mRNA changes after 2.5 h of cycling: influence of carbohydrate. *Med. Sci. Sports Exerc.* **37**, 1283–1290 (2005).
54. Bastaki, M. *et al.* Genotype-activity relationship for Mn-superoxide dismutase, glutathione peroxidase 1 and catalase in humans. *Pharmacogenet. Genomics* **16**, 279–286 (2006).
55. García-López, D. *et al.* Effects of strength and endurance training on antioxidant enzyme gene expression and activity in middle-aged men. *Scand. J. Med. Sci. Sports* **17**, 595–604 (2007).
56. Mastaloudis, A., Morrow, J. D., Hopkins, D. W., Devaraj, S. & Traber, M. G. Antioxidant supplementation prevents exercise-induced lipid peroxidation, but not inflammation, in ultramarathon runners. *Free Radic. Biol. Med.* **36**, 1329–1341 (2004).
57. Mastaloudis, A. *et al.* Endurance exercise results in DNA damage as detected by the comet assay. *Free Radic. Biol. Med.* **36**, 966–975 (2004).
58. Viktor Bauer, R. S. Bauer V, Sotnikova R. Nitric oxide—the endothelium-derived relaxing factor and its role in endothelial functions. *Gen. Physiol. Biophys.* **29**, 319–40 (2010).
59. Puddu, P., Cravero, E., Puddu, G. M. & Muscari, A. Genes and atherosclerosis: at the origin of the predisposition. *Int. J. Clin. Pract.* **59**, 462–472 (2005).
60. Ferguson, J. F. *et al.* NOS3 gene polymorphisms are associated with risk markers of cardiovascular disease, and interact with omega-3 polyunsaturated fatty acids. *Atherosclerosis* **211**, 539–544 (2010).
61. van Beynum, I. M. *et al.* Common 894G>T single nucleotide polymorphism in the gene coding for endothelial nitric oxide synthase (eNOS) and risk of congenital heart defects. *Clin. Chem. Lab. Med.* **46**, 1369–1375 (2008).

62. Hao, Y., Montiel, R. & Huang, Y. Endothelial nitric oxide synthase (eNOS) 894 G>T polymorphism is associated with breast cancer risk: a meta-analysis. *Breast Cancer Res. Treat.* **124**, 809–813 (2010).
63. Jang, M. J. *et al.* Association of eNOS polymorphisms (-786T>C, 4a4b, 894G>T) with colorectal cancer susceptibility in the Korean population. *Gene* **512**, 275–281 (2013).
64. Montgomery, H. E. *et al.* Association of angiotensin-converting enzyme gene I/D polymorphism with change in left ventricular mass in response to physical training. *Circulation* **96**, 741–747 (1997).
65. Corvol, P. & Jeunemaitre, X. Molecular Genetics of Human Hypertension: Role of Angiotensinogen. *Endocr. Rev.* **18**, 662–677 (1997).
66. Gomez-Gallego, F. *et al.* The C allele of the AGT Met235Thr polymorphism is associated with power sports performance. *Appl. Physiol. Nutr. Metab. Physiol. Appl. Nutr. Metab.* **34**, 1108–1111 (2009).
67. Zarębska, A. *et al.* Association of rs699 (M235T) polymorphism in the AGT gene with power but not endurance athlete status. *J. Strength Cond. Res.* **27**, 2898–2903 (2013).
68. Aleksandra, Z. *et al.* The AGT Gene M235T Polymorphism and Response of Power-Related Variables to Aerobic Training. *J. Sports Sci. Med.* **15**, 616–624 (2016).
69. Paillard, F. *et al.* Genotype-phenotype relationships for the renin-angiotensin-aldosterone system in a normal population. *Hypertens. Dallas Tex 1979* **34**, 423–429 (1999).
70. Rankinen, T. *et al.* AGT M235T and ACE ID polymorphisms and exercise blood pressure in the HERITAGE Family Study. *Am. J. Physiol. Heart Circ. Physiol.* **279**, H368-374 (2000).
71. Rauramaa, R. *et al.* Physical exercise and blood pressure with reference to the angiotensinogen M235T polymorphism. *Physiol. Genomics* **10**, 71–77 (2002).
72. Karjalainen, J. *et al.* Angiotensinogen gene M235T polymorphism predicts left ventricular hypertrophy in endurance athletes. *J. Am. Coll. Cardiol.* **34**, 494–499 (1999).
73. Alves, G. B. *et al.* Influence of angiotensinogen and angiotensin-converting enzyme polymorphisms on cardiac hypertrophy and improvement on maximal aerobic capacity caused by exercise training. *Eur. J. Cardiovasc. Prev. Rehabil. Off. J. Eur. Soc. Cardiol. Work. Groups Epidemiol. Prev. Card. Rehabil. Exerc. Physiol.* **16**, 487–492 (2009).
74. Diet, F. *et al.* ACE and angiotensinogen gene genotypes and left ventricular mass in athletes. *Eur. J. Clin. Invest.* **31**, 836–842 (2001).
75. Jones, A., Montgomery, H. E. & Woods, D. R. Human performance: a role for the ACE genotype? *Exerc. Sport Sci. Rev.* **30**, 184–190 (2002).
76. Williams, A. G. *et al.* Circulating angiotensin converting enzyme activity is correlated with muscle strength. *Med. Sci. Sports Exerc.* **37**, 944–948 (2005).

77. Charbonneau, D. E. *et al.* ACE genotype and the muscle hypertrophic and strength responses to strength training. *Med. Sci. Sports Exerc.* **40**, 677–683 (2008).
78. Zhang, B. *et al.* The I allele of the angiotensin-converting enzyme gene is associated with an increased percentage of slow-twitch type I fibers in human skeletal muscle. *Clin. Genet.* **63**, 139–144 (2003).
79. Myerson, S. *et al.* Human angiotensin I-converting enzyme gene and endurance performance. *J. Appl. Physiol. Bethesda Md* **1985** **87**, 1313–1316 (1999).
80. Woods, D. *et al.* Elite swimmers and the D allele of the ACE I/D polymorphism. *Hum. Genet.* **108**, 230–232 (2001).
81. Costa, A. M. *et al.* Association between ACE D allele and elite short distance swimming. *Eur. J. Appl. Physiol.* **106**, 785–790 (2009).
82. Nazarov, I. B. *et al.* The angiotensin converting enzyme I/D polymorphism in Russian athletes. *Eur. J. Hum. Genet. EJHG* **9**, 797–801 (2001).
83. Papadimitriou, I. D., Papadopoulos, C., Kouvatsi, A. & Triantaphyllidis, C. The ACE I/D polymorphism in elite Greek track and field athletes. *J. Sports Med. Phys. Fitness* **49**, 459–463 (2009).
84. Ahmetov, I. I. *et al.* The association of ACE, ACTN3 and PPARA gene variants with strength phenotypes in middle school-age children. *J. Physiol. Sci. JPS* **63**, 79–85 (2013).
85. Giaccaglia, V. *et al.* Interaction between angiotensin converting enzyme insertion/deletion genotype and exercise training on knee extensor strength in older individuals. *Int. J. Sports Med.* **29**, 40–44 (2008).
86. Pereira, A. *et al.* ACE I/D and ACTN3 R/X polymorphisms as potential factors in modulating exercise-related phenotypes in older women in response to a muscle power training stimuli. *Age Dordr. Neth.* **35**, 1949–1959 (2013).
87. Ahmetov, I. I. & Rogozkin, V. A. Genes, athlete status and training -- An overview. *Med. Sport Sci.* **54**, 43–71 (2009).
88. Gayagay, G. *et al.* Elite endurance athletes and the ACE I allele--the role of genes in athletic performance. *Hum. Genet.* **103**, 48–50 (1998).
89. Cieszczyk, P., Krupecki, K., Maciejewska, A. & Sawczuk, M. The angiotensin converting enzyme gene I/D polymorphism in Polish rowers. *Int. J. Sports Med.* **30**, 624–627 (2009).
90. Alvarez, R. *et al.* Genetic variation in the renin-angiotensin system and athletic performance. *Eur. J. Appl. Physiol.* **82**, 117–120 (2000).
91. Collins, M. *et al.* The ACE gene and endurance performance during the South African Ironman Triathlons. *Med. Sci. Sports Exerc.* **36**, 1314–1320 (2004).
92. Shenoy, S., Tandon, S., Sandhu, J. & Bhanwer, A. S. Association of Angiotensin Converting Enzyme gene Polymorphism and Indian Army Triathletes Performance. *Asian J. Sports Med.* **1**, 143–150 (2010).
93. Lucía, A. *et al.* Is there an association between ACE and CKMM polymorphisms and cycling performance status during 3-week races? *Int. J. Sports Med.* **26**, 442–447 (2005).

94. Hruskovicová, H. *et al.* The angiotensin converting enzyme I/D polymorphism in long distance runners. *J. Sports Med. Phys. Fitness* **46**, 509–513 (2006).
95. Min, S.-K. *et al.* Is there a gender difference between ACE gene and race distance? *Appl. Physiol. Nutr. Metab. Physiol. Appl. Nutr. Metab.* **34**, 926–932 (2009).
96. Znazen, H. *et al.* Genetic advantageous predisposition of angiotensin converting enzyme id polymorphism in Tunisian athletes. *J. Sports Med. Phys. Fitness* **56**, 724–730 (2016).
97. Ma, F. *et al.* The association of sport performance with ACE and ACTN3 genetic polymorphisms: a systematic review and meta-analysis. *PLoS One* **8**, e54685 (2013).
98. Papadimitriou, I. D. *et al.* ACTN3 R577X and ACE I/D gene variants influence performance in elite sprinters: a multi-cohort study. *BMC Genomics* **17**, (2016).
99. Williams, A. G. *et al.* Bradykinin receptor gene variant and human physical performance. *J. Appl. Physiol. Bethesda Md* **1985** **96**, 938–942 (2004).
100. Saunders, C. J. *et al.* The bradykinin beta 2 receptor (BDKRB2) and endothelial nitric oxide synthase 3 (NOS3) genes and endurance performance during Ironman Triathlons. *Hum. Mol. Genet.* **15**, 979–987 (2006).
101. Tsianos, G. I. *et al.* Associations of polymorphisms of eight muscle- or metabolism-related genes with performance in Mount Olympus marathon runners. *J. Appl. Physiol. Bethesda Md* **1985** **108**, 567–574 (2010).
102. Saunders, C. J. *et al.* Dipsogenic genes associated with weight changes during Ironman Triathlons. *Hum. Mol. Genet.* **15**, 2980–2987 (2006).
103. Sawczuk, M. *et al.* The -9/+9 polymorphism of the bradykinin receptor Beta 2 gene and athlete status: a study involving two European cohorts. *Hum. Biol.* **85**, 741–756 (2013).
104. Grenda, A., Leońska-Duniec, A., Cięższyk, P. & Zmijewski, P. Bdkrb2 gene -9/+9 polymorphism and swimming performance. *Biol. Sport* **31**, 109–113 (2014).
105. Prior, S. J. *et al.* DNA sequence variation in the promoter region of the VEGF gene impacts VEGF gene expression and maximal oxygen consumption. *Am. J. Physiol. Heart Circ. Physiol.* **290**, H1848-1855 (2006).
106. Booth, F. W., Rueggsegger, G. N., Toedebusch, R. G. & Yan, Z. Endurance Exercise and the Regulation of Skeletal Muscle Metabolism. *Prog. Mol. Biol. Transl. Sci.* **135**, 129–151 (2015).
107. Coffey, V. G. *et al.* Interaction of contractile activity and training history on mRNA abundance in skeletal muscle from trained athletes. *Am. J. Physiol. Endocrinol. Metab.* **290**, E849-855 (2006).
108. Gustafsson, T., Puntchart, A., Kaijser, L., Jansson, E. & Sundberg, C. J. Exercise-induced expression of angiogenesis-related transcription and growth factors in human skeletal muscle. *Am. J. Physiol.* **276**, H679-685 (1999).

109. Watson, C. J., Webb, N. J., Bottomley, M. J. & Brenchley, P. E. Identification of polymorphisms within the vascular endothelial growth factor (VEGF) gene: correlation with variation in VEGF protein production. *Cytokine* **12**, 1232–1235 (2000).
110. Akhmetov, I. I. *et al.* [Polymorphism of the vascular endothelial growth factor gene (VEGF) and aerobic performance in athletes]. *Fiziol. Cheloveka* **34**, 97–101 (2008).
111. Ahmetov, I. I. *et al.* The combined impact of metabolic gene polymorphisms on elite endurance athlete status and related phenotypes. *Hum. Genet.* **126**, 751–761 (2009).
112. He, Z. *et al.* NRF2 genotype improves endurance capacity in response to training. *Int. J. Sports Med.* **28**, 717–721 (2007).
113. Eynon, N. *et al.* NRF2 intron 3 A/G polymorphism is associated with endurance athletes' status. *J. Appl. Physiol. Bethesda Md 1985* **107**, 76–79 (2009).
114. Maciejewska-Karłowska, A. *et al.* The GABPB1 gene A/G polymorphism in Polish rowers. *J. Hum. Kinet.* **31**, 115–120 (2012).
115. Ling, C. *et al.* Multiple environmental and genetic factors influence skeletal muscle PGC-1alpha and PGC-1beta gene expression in twins. *J. Clin. Invest.* **114**, 1518–1526 (2004).
116. Norrbom, J. *et al.* PGC-1alpha mRNA expression is influenced by metabolic perturbation in exercising human skeletal muscle. *J. Appl. Physiol. Bethesda Md 1985* **96**, 189–194 (2004).
117. Lucia, A. *et al.* PPARGC1A genotype (Gly482Ser) predicts exceptional endurance capacity in European men. *J. Appl. Physiol. Bethesda Md 1985* **99**, 344–348 (2005).
118. Stefan, N. *et al.* Genetic variations in PPARD and PPARGC1A determine mitochondrial function and change in aerobic physical fitness and insulin sensitivity during lifestyle intervention. *J. Clin. Endocrinol. Metab.* **92**, 1827–1833 (2007).
119. Akhmetov, I. I. *et al.* [Association of regulatory genes polymorphisms with aerobic and anaerobic performance of athletes]. *Ross. Fiziol. Zh. Im. I M Sechenova* **93**, 837–843 (2007).
120. Eynon, N. *et al.* Do PPARGC1A and PPARalpha polymorphisms influence sprint or endurance phenotypes? *Scand. J. Med. Sci. Sports* **20**, e145-150 (2010).
121. Maciejewska, A., Sawczuk, M., Cieszczyk, P., Mozhayskaya, I. A. & Ahmetov, I. I. The PPARGC1A gene Gly482Ser in Polish and Russian athletes. *J. Sports Sci.* **30**, 101–113 (2012).
122. Gineviciene, V. *et al.* Association analysis of ACE, ACTN3 and PPARGC1A gene polymorphisms in two cohorts of European strength and power athletes. *Biol. Sport* **33**, 199–206 (2016).
123. Ahmetov, I. I. *et al.* PPARalpha gene variation and physical performance in Russian athletes. *Eur. J. Appl. Physiol.* **97**, 103–108 (2006).

124. Jamshidi, Y. *et al.* Peroxisome proliferator-activated receptor alpha gene regulates left ventricular growth in response to exercise and hypertension. *Circulation* **105**, 950–955 (2002).
125. Kochanska-Dziurawicz, A. A. *et al.* CATECHOLAMINES AND β 2-ADRENOCEPTOR GENE EXPRESSION BEFORE AND AFTER MAXIMAL INCREMENTAL CYCLE TEST IN YOUNG ICE HOCKEY PLAYERS: RELATION TO WORK PERFORMED. *Biol. Sport* **30**, 85–90 (2013).
126. Snyder, E. M. *et al.* Arg16Gly polymorphism of the beta2-adrenergic receptor is associated with differences in cardiovascular function at rest and during exercise in humans. *J. Physiol.* **571**, 121–130 (2006).
127. Eisenach, J. H., McGuire, A. M., Schwingler, R. M., Turner, S. T. & Joyner, M. J. The Arg16/Gly beta2-adrenergic receptor polymorphism is associated with altered cardiovascular responses to isometric exercise. *Physiol. Genomics* **16**, 323–328 (2004).
128. Wolfarth, B. *et al.* Association between a beta2-adrenergic receptor polymorphism and elite endurance performance. *Metabolism.* **56**, 1649–1651 (2007).
129. Wagoner, L. E. *et al.* Polymorphisms of the beta(2)-adrenergic receptor determine exercise capacity in patients with heart failure. *Circ. Res.* **86**, 834–840 (2000).
130. Moore, G. E. *et al.* Obesity gene variant and elite endurance performance. *Metabolism.* **50**, 1391–1392 (2001).
131. Sawczuk, M., Maciejewska-Karlowska, A., Cieszczyk, P., Skotarczak, B. & Ficek, K. Association of the ADRB2 Gly16Arg and Glu27Gln polymorphisms with athlete status. *J. Sports Sci.* **31**, 1535–1544 (2013).
132. Macho-Azcarate, T., Marti, A., Calabuig, J. & Martinez, J. A. Basal fat oxidation and after a peak oxygen consumption test in obese women with a beta2 adrenoceptor gene polymorphism. *J. Nutr. Biochem.* **14**, 275–279 (2003).
133. Liu, X.-G. *et al.* Genome-wide Association and Replication Studies Identified TRHR as an Important Gene for Lean Body Mass. *Am. J. Hum. Genet.* **84**, 418–423 (2009).
134. Yang, N. *et al.* ACTN3 genotype is associated with human elite athletic performance. *Am. J. Hum. Genet.* **73**, 627–631 (2003).
135. Eynon, N. *et al.* Is the interaction between HIF1A P582S and ACTN3 R577X determinant for power/sprint performance? *Metabolism.* **59**, 861–865 (2010).
136. Roth, S. M. *et al.* The ACTN3 R577X nonsense allele is under-represented in elite-level strength athletes. *Eur. J. Hum. Genet. EJHG* **16**, 391–394 (2008).
137. Gómez-Gallego, F. *et al.* Endurance performance: genes or gene combinations? *Int. J. Sports Med.* **30**, 66–72 (2009).
138. Shang, X., Huang, C., Chang, Q., Zhang, L. & Huang, T. Association between the ACTN3 R577X polymorphism and female endurance athletes in China. *Int. J. Sports Med.* **31**, 913–916 (2010).
139. Döring, F. E. *et al.* ACTN3 R577X and other polymorphisms are not associated with elite endurance athlete status in the Genathlete study. *J. Sports Sci.* **28**, 1355–1359 (2010).

140. Saunders, C. J. *et al.* No association of the ACTN3 gene R577X polymorphism with endurance performance in Ironman Triathlons. *Ann. Hum. Genet.* **71**, 777–781 (2007).
141. Lucia, A. *et al.* ACTN3 genotype in professional endurance cyclists. *Int. J. Sports Med.* **27**, 880–884 (2006).
142. Ogan, D. & Pritchett, K. Vitamin D and the athlete: risks, recommendations, and benefits. *Nutrients* **5**, 1856–1868 (2013).
143. Geusens, P. *et al.* Quadriceps and grip strength are related to vitamin D receptor genotype in elderly nonobese women. *J. Bone Miner. Res. Off. J. Am. Soc. Bone Miner. Res.* **12**, 2082–2088 (1997).
144. Grundberg, E. *et al.* Genetic variation in the human vitamin D receptor is associated with muscle strength, fat mass and body weight in Swedish women. *Eur. J. Endocrinol. Eur. Fed. Endocr. Soc.* **150**, 323–328 (2004).
145. Wang, P. *et al.* Association between polymorphisms of vitamin D receptor gene Apal, Bsm1 and Taq1 and muscular strength in young Chinese women. *Int. J. Sports Med.* **27**, 182–186 (2006).
146. Windelinckx, A. *et al.* Polymorphisms in the vitamin D receptor gene are associated with muscle strength in men and women. *Osteoporos. Int. J. Establ. Result Coop. Eur. Found. Osteoporos. Natl. Osteoporos. Found. USA* **18**, 1235–1242 (2007).
147. Bahat, G. *et al.* Bsm1 polymorphism in the vitamin D receptor gene is associated with leg extensor muscle strength in elderly men. *Aging Clin. Exp. Res.* **22**, 198–205 (2010).
148. Ceglia, L. & Harris, S. S. Vitamin D and its role in skeletal muscle. *Calcif. Tissue Int.* **92**, 151–162 (2013).
149. Bao, L. *et al.* Association between vitamin D receptor Bsm1 polymorphism and bone mineral density in pediatric patients: A meta-analysis and systematic review of observational studies. *Medicine (Baltimore)* **96**, e6718 (2017).
150. Womack, C. J. *et al.* The influence of a CYP1A2 polymorphism on the ergogenic effects of caffeine. *J. Int. Soc. Sports Nutr.* **9**, 7 (2012).
151. Cornelis, M. C., El-Sohemy, A., Kabagambe, E. K. & Campos, H. Coffee, CYP1A2 genotype, and risk of myocardial infarction. *JAMA* **295**, 1135–1141 (2006).
152. Pataký, M. W. *et al.* Caffeine and 3-km cycling performance: Effects of mouth rinsing, genotype, and time of day. *Scand. J. Med. Sci. Sports* **26**, 613–619 (2016).
153. Salinero, J. J. *et al.* CYP1A2 Genotype Variations Do Not Modify the Benefits and Drawbacks of Caffeine during Exercise: A Pilot Study. *Nutrients* **9**, (2017).
154. Goel, N., Basner, M., Rao, H. & Dinges, D. F. Circadian Rhythms, Sleep Deprivation, and Human Performance. *Prog. Mol. Biol. Transl. Sci.* **119**, 155–190 (2013).
155. Vitale, J. A. & Weydahl, A. Chronotype, Physical Activity, and Sport Performance: A Systematic Review. *Sports Med. Auckl. NZ* (2017). doi:10.1007/s40279-017-0741-z
156. Katzenberg, D. *et al.* A CLOCK polymorphism associated with human diurnal preference. *Sleep* **21**, 569–576 (1998).

157. Garaulet, M. *et al.* Ghrelin, Sleep Reduction and Evening Preference: Relationships to CLOCK 3111 T/C SNP and Weight Loss. *PLoS ONE* **6**, (2011).
158. Mishima, K., Tozawa, T., Satoh, K., Saitoh, H. & Mishima, Y. The 3111T/C polymorphism of hClock is associated with evening preference and delayed sleep timing in a Japanese population sample. *Am. J. Med. Genet. Part B Neuropsychiatr. Genet. Off. Publ. Int. Soc. Psychiatr. Genet.* **133B**, 101–104 (2005).
159. Rae, D. E., Stephenson, K. J. & Roden, L. C. Factors to consider when assessing diurnal variation in sports performance: the influence of chronotype and habitual training time-of-day. *Eur. J. Appl. Physiol.* **115**, 1339–1349 (2015).